

Trade Terms	FOB Shenzhen / Guangzhou / Hongkong
Production Capacity	500,000PCS/Month
Payment	T/T, PayPal.
Delivery Time	3-7 working days for making sample. 10 to 15 working days for production.
Packing	Standard Export Cartons/Custom Cartons
Shipping	By sea, By air or by express delivery (DHL/UPS/FEDEX)
MOQ	100PCS
Place of Origin	Dongguan, Guangzhou, China
Certification	CE,FCC,ROHS,REACH,UKCA

Product Features:

If your doctor has recommended transcutaneous vagus nerve stimulation (tVNS) you will need a good fitting ear electrode.

The tragus clip electrode has been designed by experienced health professionals to be a comfortable, adjustable, and secure fit.

We are currently working through the relevant regulatory approvals processes,

We make many different kinds of CES ear clips electrodes, One sided with conductive rubber ear clip electrodes, Double sided ear clip electrodes and also wired CES Ear clip electrodes so customer can use it in different machines and conditions, If you need to customize one Tragus Electrodes ear clip please feel free to contact with us we will give you professional advise and design,



These Ear Clip Electrodes are a must-have for effective Cranial Electrotherapy Stimulation (CES), by applying to different points of the human ear, along with the use of a TENS or Microcurrent therapy device, to treat such conditions as:

Anxiety,
Depression,
Post-Traumatic Stress Disorder (PTSD),
Mood Swings,
TMJ,
Tinnitus,
Insomnia,
Acute & Chronic Pain,
Other Neurological & Psychiatric Conditions.

Ear Clip Electrodes are commonly placed on the ear lobe, but this specific style is also optimal for placement on the tragus or cymba conchae of the left ear for purposes of Transcutaneous Vagus

Nerve Stimulation (tVNS).

Features one 1/4" diameter conductive rubber electrode on one end of the Clip, which is recommended to be used with a conductive gel, spray, or water, before applying the Ear Clip, to ensure the most effective and comfortable stimulation.

Very sturdy and stable spring design to help keep the Clip in place on the ear, compared to other similar items in the market.

Inside of the Clip itself is an industry-standard, 2mm-diameter female pin-style connection for simple connectivity and compatibility with most electrotherapy lead wires and devices.





Company Introduction

[Goochain Technology Co.,Ltd](#) is a professional custom cable assemblies manufacturer, manufacturing custom usb cables, Magnetic pogo pin cables manufacturer , Medical ecg cables ,targus clip electrodes supplier, PoweredUSB extension Cables, Ethernet cables,Medical plugs and chargers,Wiring harnesses and computer cables , Founded in year of 2015 and we have 2000 square meters custom cable assemblies workshop, around 80 workers, 1 shift is 8 hours per day. We followed international ISO management Standard to make our custom cables, medical cables and chargers with Chinese factory price .



Certificates:

Goochain custom usb cables, Medical cables, computer cables, Pogo pin cables, USB Type c cables, USB Type C connectors, PD Fast chargers, Powered USB extension Cables all has passed the international required certificate so we can export all over the world,



FAQ

Q1: Can I start a sample order first before ordering more?

A: Sure, Goochain provide free samples if we have stock, sample customized should be charged some set up cost. We appreciate you to order sample from us to check our quality.

Q2: What payment methods do you accept?

A: Alibaba Trade Assurance, PayPal, T/T, Bank Transfer, Western Union and L/C

Q3: How does your factory do regarding quality control?

A: Quality is priority. QC department with professional quality assurance specialist attach great importance to quality controlling from the very beginning to the very end. Our factory has gained Rohs, SGS, CE FCC,ISO9001:2008 authentication. Etc

Q4: Factory? Yes! Can I customize items I like, such as to print our own logo, or to alter the design?

A: Yes, we are a factory with more than 100 employees. We offer OEM/ODM service. Please contact us or our sales representative to check with the details.

Q5: Delivery time and warranty? Which shipping service you will use for my goods?

A: Our lead time is around 2-4 weeks, and we can provide 12 months warranty since delivery. We used DHL,Fedex,UPS,TNT,EMS as customer like.

Q6: Where is your factory located? How can I visit there?

A: Sure. We are more than happy to receive you in Dongguan and show you around our workplace. Our factory is located in 3rd. Floor,Building A, No.8 Kangyang 1st Street,Shebei Village,Huangjiang Town,Dongguan City,Guangdong, China.